

Latex Cover	Model #	Arm/Leg Circumference	Length
X-Small Full Arm	<b>FA-12</b>	6 – 6.75 inches	16 inches
Small Full Arm	<b>FA-14</b>	6.75 – 8.75 inches	23 inches
Small Half Arm	<b>HA-13</b>	7.75 – 10 inches	17.5 inches
Large Half Arm	<b>HA-15</b>	10 inches & up	20 inches
Medium Full Arm	<b>FA-16</b>	8.75 – 10 inches	28 inches
Large Full Arm	<b>FA-18</b>	10 inches & up	31.5 inches
<b>X-Small Full Leg</b>			
X-Small Full Leg	<b>FL-12</b>	7.5 – 11 inches	19 inches
Small Half Leg	<b>HL-13</b>	10 – 13 inches	21 inches
Large Half Leg	<b>HL-15</b>	13 inches & up	23.5 inches
Small Full Leg	<b>FL-14</b>	14 – 16.5 inches	29 inches
Medium Full Leg	<b>FL-16</b>	16.5 – 21 inches	33 inches
Large Full Leg	<b>FL-18</b>	21 inches & up	37 inches
<b>*Keep in mind, to achieve the best fit, consider both half and full sizes.</b>			
<b>Measuring Circumference:</b> Measure the circumference of the arm or leg at 2 inches above the cast. If the end of the cast is very close to the elbow/knee, then measure mid-bicep/mid-thigh. Do not measure around the cast! Only measure around the skin of the arm/leg above the cast.			
<b>Measuring Length:</b> Measure from the tip of the middle finger or heel of the foot to the top end of the cast. It is okay if the cast protector measurement is longer than the limb measurement. The protector will still work. Extra length will be removed when the vacuum seal is created.			
<b>X-Small Prosthetic Leg</b>			
X-Small Prosthetic Leg	<b>PL-11</b>	8 -12 inches	29 inches
Small Prosthetic Leg	<b>PL-13</b>	12 – 16 inches	32 inches
Medium Prosthetic Leg	<b>PL-15</b>	15 – 19 inches	35 inches
Large Prosthetic Leg	<b>PL-17</b>	17 – 21 inches	38 inches
X-Large Prosthetic Leg	<b>PL-19</b>	21 inches & Up	41 inches
<b>Stretch on the DryPro Prosthetic Leg Cover:</b> Pull over the leg to see how it fits. If the DryPro is too tight, you can trim 1/4 inch along the lines towards the top of the opening and try it on again. Be sure not to trim too much, as a tight seal is needed for the vacuum to work. Stretch the DryPro over the prosthetic so the top fits smooth against the skin. Allow some slack as the vacuum will remove the excess material.			

**Call us at 1-888-337-9776**

[www.DRYCORP.com](http://www.DRYCORP.com)